
Lab

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
Lab

Laboport® N96
TRANSLATION OF ORIGINAL
OPERATING INSTRUCTION
ENGLISH

EN

LABOPORT®

LABORATORY PUMP

Notice!

Before operating the pump and the accessories, read the operating and installation instructions and observe the safety notices!

KNF develops, produces and sells high-quality diaphragm pumps and systems for gases, vapors and fluids. They are available as built-in pumps (OEM), for process applications and for use in laboratories.

With an extensive range application-specific pumps and systems, KNF has been a technology leader for decades. Countless customers around the world trust KNF as a supplier of long-lasting and reliable products for a wide array of different applications.

Thank you for placing your trust in us and for deciding to purchase one of our products.

Quality has top priority with our products, which is why each of our products undergoes a 100% test. All processes are certified in accordance with ISO 9001, and the RoHS- and Reach guidelines are satisfied for all materials.

KNF Neuberger GmbH
Alter Weg 3
79112 Freiburg
Germany
Tel. +49 (0)7664/5909-0

E-mail: info.de@knf.com
www.knf.com

Index

- 1 Scope of delivery 4
- 2 About this document 5
 - 2.1 Using the operating instructions..... 5
 - 2.2 Exclusion of liability..... 5
 - 2.3 Applicable documents..... 6
 - 2.4 Symbols and markings..... 6
- 3 Use 9
 - 3.1 Proper use 9
 - 3.2 Improper use 10
- 4 Safety..... 11
- 5 Technical data 14
- 6 Design and function 16
- 7 Transport 18
 - General 18
- 8 Set up, installation and connection 20
 - 8.1 Connecting the pump 20
- 9 Operation 23
 - 9.1 Pump..... 23
 - 9.2 Information on switching the pump on and off
..... 26
- 10 Servicing 28
 - 10.1 Servicing schedule..... 28
 - 10.2 Cleaning..... 29
 - 10.3 Replacing diaphragm and valve plates 31
- 11 Spare parts and accessories 40
 - 11.1 Spare parts 40
 - 11.2 Accessories..... 41
- 12 Troubleshooting 42
- 13 Returns 45
- 14 Index 46

1 Scope of delivery

- Laboport® N96
- Power supply incl. plug insert (EU, US, UK, AU)
- Hose connector with NPT thread ID6
- Silencer
- Operating instructions


Fig. 1: Scope of delivery

Unpacking the pump

Inspect the pump and the included accessories for transport damage after unpacking. If the packaging is damaged, inform the responsible forwarder so that a damage report can be prepared. For further information, read Chapter 7 *Transport*.

2 About this document

2.1 Using the operating instructions

The operating instructions are part of the pump.

- In the event of uncertainties with regard to the content of the operating instructions, please contact the manufacturer (contact data: see www.knf.com). Please have the type designation and serial number of the pump ready.
- Read the operating instructions before you commission the pump.
- Give the operating instructions to the next owner.
- Keep the operating instructions within reach at all times.

Project pumps For customer-specific project pumps (pump models that begin with "PJ" or "PM"), there may be deviations from the operating and installation instructions.

- For project pumps, also observe the agreed specifications.

2.2 Exclusion of liability

The manufacturer assumes no liability for damages and malfunctions resulting from failure to observe the operating instructions.

The manufacturer assumes no liability for damages and malfunctions resulting from changes or modifications to the device and improper handling.

The manufacturer assumes no liability for damages and malfunctions resulting from impermissible spare parts and accessories.

2.3 Applicable documents

The listed documents must also be observed. The valid versions are available at www.knf.com/Downloads.

- Data sheet
- 3D model

Also to be observed are:

- Local terms and conditions
- Sales documents and agreement between KNF and customer.

2.4 Symbols and markings

Warning notice

A notice that warns you of danger is located here.

Possible consequences of a failure to observe the warning notice are specified here. The signal word, e.g., warning, indicates the danger level.

- Measures for avoiding the danger and its consequences are specified here.
-

Danger levels

Signal word	Meaning	Consequences if not observed
DANGER	warns of immediate danger	Death or serious injury or serious damage will result.
WARNING	warns of possible danger	Death or serious injury or serious damage is possible.
CAUTION	warns of a possibly dangerous situation	Minor injuries or damage is possible.
NOTICE	Warns of possible damage	Damage is possible.

*Tab. 1: Danger levels***Other notices and symbols**

→ An activity to be carried out is specified here (a step).

1. The first step of an activity to be carried out is specified here.
Follow other sequentially numbered steps.

i This symbol indicates important information.

Explanations of pictograms

Pictogram	Meaning
	General warning symbol
	Warning of hot surface
	Warning of electrical voltage
	Warning of explosive materials
	Warning of poisonous substances
	Observe the operating instructions
	General mandatory sign
	Unplug mains plug
	Use foot protection
	Use hand protection
	Environmentally conscious disposal
	Recycling

Tab.2: Explanations of pictograms

3 Use

3.1 Proper use

The pumps are intended exclusively for transferring gases and vapors.

Responsibility of the owner

Operating parameters and conditions

Only install and operate the pumps in accordance with the operating parameters and conditions described in Chapter 5 *Technical data*.

Only pumps that are fully assembled and in the condition as delivered may be operated.

Make sure that the installation location is dry and that the pump is protected against rain, splash, gushing, and drip water as well as from other contaminants.

Check the tightness of the connections between the pipes of the application and the pump (or the connection of the pump) at regular intervals. Leaky connections carry the risk of releasing dangerous gases and vapors from the pump system.

Requirements on the transferred medium

Before transferring a medium, check whether the medium can be transferred danger-free in the specific application.

Before using a medium, check the compatibility of the media-contacting components (see 5 *Technical data*) with the medium.

Risk of dangerous gas mixtures during pump operation if diaphragm breaks: Depending on the medium being transferred, breakage of the diaphragm can result in a dangerous mixture if the medium mixes with the air in the compressor housing or the surroundings.

Only transfer gases that remain stable under the pressures and temperatures that arise in the pump.

Accessories

Laboratory accessories or additional components that are connected to a pump must be designed for the pneumatic data of the pump (see 5 *Technical data*).

3.2 Improper use

The pumps are not allowed to be operated in explosive atmospheres.

The pumps are not suitable for delivering:

- Dusts
- Fluids
- Aerosols
- Biological and microbiological substances
- Fuel
- Explosive substances and flammable materials
- Fibers
- Oxidants
- Food

Pumps that can generate both a vacuum and overpressure must not be used for the simultaneous generation of a vacuum and overpressure.

This function can be made possible on a project basis following consultation with KNF Customer Service.

No overpressure may be applied to the suction side of the pump.

This function can be made possible on a project basis following consultation with KNF Customer Service.

4 Safety

i Observe the safety notices in Chapters Installation and connection and Operation.

The pumps are built in accordance with the generally recognized rules of engineering and the occupational health and safety and accident prevention regulations. Nevertheless, dangers can arise during their use that lead to injuries to the user or third parties or to damage to the pump or other property.

Only use the pumps in perfect technical condition, for their intended purpose, safely and aware of the dangers and in observation of the operating and installation instructions.

The components that are to be connected to the pumps must be designed according to the pneumatic data of the pumps.

When connecting the pumps to the electrical mains, observe the corresponding safety rules.

Personnel Make sure that only specially trained and instructed personnel or specialist personnel work on the pumps. This applies, in particular, to assembly, connection and servicing work.

Make sure that the personnel have read and understood the operating and installation instructions, particularly the chapter on safety.

Working in a safety conscious manner Observe the regulations on accident prevention and safety during all work on the pumps and during operation.

Avoid contact with the heads and housing parts, as the pump heats up during operation.

When working on the pump, make certain that the pump is disconnected from mains and without power.

Make certain that no dangers arise from flows when gas connections are open, from the effects of noise or from hot, corrosive, dangerous and environmentally hazardous gases.

Make sure that an EMC-compliant installation of the pump is ensured at all times and that no dangerous situation can thereby arise.

Working with hazardous media When pumping hazardous media, observe the safety regulations for the handling of these media.

**Working with
combustible
media**

Note that the pumps are not designed to be explosion-proof.

Make certain that the temperature of the medium is always sufficiently below the ignition temperature of the medium so as to prevent ignition or explosion. This also applies for abnormal operating situations.

Note here that the temperature of the medium increases when the pump compresses the medium.

Therefore, make certain that the temperature of the medium also remains sufficiently below the ignition temperature of the medium even when it is compressed to the maximum permissible operating pressure of the pump. The maximum permissible operating pressure of the pump is stated in Chapter 5 *Technical data*.

Make certain that the permissible ambient temperature (5 *Technical data*) is not exceeded.

If applicable, also take into consideration external energy sources (e.g., radiation sources) that could additionally heat the medium.

In case of doubt, contact KNF Customer Service.

**Environmental
protection**

Store and dispose of the pump as well as all replacement parts in accordance with environmental regulations. Observe the respective national and international regulations. This applies in particular to parts that are contaminated with toxic substances.

Dispose of packaging material that is no longer needed in an environmentally responsible manner. The packaging materials are recyclable.

Dispose of old devices in an environmentally responsible manner. Dispose of old devices via suitable collection systems. Old devices contain valuable recyclable materials.

EU/EC directives/standards

- The pumps correspond to directives:
- 2011/65/EU (RoHS)
 - 2014/30/EU (EMC)
 - 2006/42/EC (MRL).

The following harmonized standards are satisfied:

- DIN EN 50581
- DIN EN 1012-2
- DIN EN ISO 12100
- DIN EN 61010-1
- DIN EN 61326-1

The pumps are inspected by TÜV according to:

- UL 61010-1
- CAN/CSA-C22.2 No. 61010-1

Certificate: No. U8 045055 0021 Rev. 00

Customer service and repairs

The pumps are maintenance-free. KNF does, however, recommend periodically inspecting the pumps for conspicuous changes in the noises and vibrations.

Only have repairs to the pumps performed by the responsible KNF Customer Service.

Housings with live components may only be opened by specialist personnel.

Use only original parts from KNF during servicing work.

5 Technical data

Technical data

Pump materials

Assembly	Material
Pump head	PPS
Diaphragm	PTFE-coated
Valves	FPM

Tab.3: Pump materials

Pneumatic data

Parameter	Value
Max. permissible operating pressure [bar rel*]	2,5
Ultimate vacuum [mbar abs.]	<130
Flow rate at atm. pressure [l/min]**	1,5 -7,0 ± 10%

Tab.4: Pneumatic data

*Bar rel related to 1013 hPa

**Liters in standard state (1013 hPa, 20°C)

Pneumatic connections

Pump type	Value
Laboport N96	NPT 1/8

Tab.5: Pneumatic connections

Electrical data

Pump

Parameter	Value
Voltage [V]	24 DC
Power P [W]	19
Max. permissible mains voltage fluctuations	± 10%
Current consumption [A]	0,9

Tab.6: Electrical data of pump

Power supply

Parameter	Value
Voltage [V]	100-240 AC
Frequency [Hz]	50/60
Current consumption [A]	0.7 A / 100 V AC 0.4 A / 230 V AC

Tab.7: Electrical data of power supply

Weight

Pump type	Value [kg]
Laboport N96	1,3

Tab.8: Weight

Other parameters

Parameter	Value
Permissible ambient temperature [°C]	+ 5 to + 40
Permissible media temperature [°C]	+ 5 to + 40
Dimensions [L x H x W] [mm]	156 x 119 x 75
Highest permissible relative air humidity of the environment	80% for temperatures to 31°C, decreasing linearly to 50% at 40°C (non-condensing).
Maximum installation altitude [m above sea level]	2000
Protection class of pump	IP40
Start up against	
- Vacuum [mbar abs.]	<130
- Pressure [bar rel*]	2,5

Tab.9:

*Bar rel related to 1013 hPa

6 Design and function

- 1 Outlet
- 2 Inlet
- 3 Rotary/push knob
- 4 Silencer
- 5 Hose connector


Fig.2: Laboport N96

The pump can be switched on and off with the rotary/push knob(3). This is also used to control the flow rate (see 9.2 *Information on switching the pump on and off*).

The pump has a vibration-free mount, thereby allowing the inlet (2) and outlet (1) to be easily moved with respect to the remainder of the housing.

Function of a diaphragm pump

- 1 Outlet valve
- 2 Inlet valve
- 3 Transfer chamber
- 4 Diaphragm
- 5 Eccentric
- 6 Connecting rod
- 7 Pump drive


Fig.3: Function of a diaphragm pump

Diaphragm pumps transfer, compress (depending on the version) and evacuate gases and vapors.

The elastic diaphragm (4) is moved up and down by the eccentric (5) and the connecting rod (6). In the downwards stroke, it aspirates the gas to be transferred via the inlet valve (2). In the upwards stroke, the diaphragm presses the medium out of the pump head via the outlet valve (1). The transfer chamber (3) is separated from the pump drive (7) by the diaphragm.

7 Transport

General

CAUTION

Personal injury and/or property damage due to incorrect or improper transport of the pump

In the event of incorrect or improper transport, the pump can fall down, be damaged or injure persons.

- Use suitable auxiliary means if necessary (carrying strap, lifting gear, etc.).
- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

CAUTION

Risk of injury from sharp edges on the packaging

There is a risk of injury from cutting on the sharp edges when grabbing corners or when opening the packaging.

- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

- Transport the pump in the original packaging to the installation location.
- Store the original packaging of the pump (e.g., for later storage).
- Inspect the pump for transport damage after receiving it.
- Document any transport damage in writing.

→ Remove any transport safeguards on the pump prior to commissioning.

Parameter

Parameter	Value
Storage temperature [°C]	+ 5 to + 40
Transport temperature [°C]	- 10 to + 60
Permissible humidity (non-condensing) [%]	30 to 85

Tab.10: Transport parameters

NOTICE

Prior to commissioning, make sure that the pump has reached the ambient temperature (5 *Technical data*).

8 Set up, installation and connection

Only connect the pumps in accordance with the operating parameters and conditions described in Chapter 5 *Technical data*.

Observe the safety notices (see Chapter Safety).

- Before connecting, store the pump at the installation location to allow it to reach the room temperature (no condensation may form).

Cooling air supply

WARNING

Danger of burning on hot surfaces

Hot surfaces could occur if the pump overheats.

- When installing the pump, make sure that sufficient cooling air in-feed and discharge is ensured.

Installation location

- Make sure that the installation location is dry and that the pump is protected against rain, splash, gushing, and drip water as well as from other contaminants.
- Select a secure location (flat surface) for the pump.
- Protect the pump from dust.
- Protect the pump from vibration, impact and external damage.

8.1 Connecting the pump

1. Remove the protective plugs from the pneumatic connections of the pump.
2. Mount suction filter or silencer accessory parts (if present).

Vacuum operation For vacuum operation, connect the hose connector (see Fig. 2/5) to the suction side and the silencer (see Fig. 2/4) or your system components to the pressure side.

Pressure operation For pressure operation, connect the hose connector to the pressure side and the silencer or your system components to the suction side.

**NOTICE**

Self-sealing NPT threads are located on the pump. When tightening the accessories, do not tighten them all the way in order to avoid damaging the threads and the pump.

-
3. Connect the suction line and the pressure line.

Connected components

→ Only connect components to the pump that are designed for the pneumatic data of the pump (see Chapter 5 *Technical data*).

Pump discharge

→ When using as a vacuum pump: Safely drain the pump discharge at the pneumatic outlet of the pump.

4. Lay the suction line and the pressure line at a downward angle to prevent condensate from running into the pump.

**NOTICE**

Hose radii selected too small can affect the stability of the pump.

**NOTICE**

Secure the pressure-side connections with a fastener (e.g., hose/pipe clamp) to prevent the hoses from slipping down from the connection.

-
5. Plug the appropriate plug insert into the power supply.

**NOTICE**

Only operate the pump with the supplied SELV power supply.

6. Connect the pump to the power supply.


Fig. 4: Connecting the plug to the pump

7. Plug the plug of the power supply into a properly installed, grounded socket.

9 Operation

9.1 Pump

9.1.1 Preparing for commissioning

Before switching on the pump, ensure the following points:

	Necessary operating requirements
Pump	→ All hoses correctly connected
Pump	→ Data of the voltage supply system are consistent with the details on the type plate of the power supply. → Pump outlet not closed or restricted.

Tab. 11: Operating requirements for commissioning

9.1.2 Commissioning

Risk of burns from hot pump parts and/or hot medium

Some pump parts may be hot during or after operation of the pump.

- Allow the pump to cool after operation.
- Take protective measures to protect against touching hot parts.

**CAUTION**

Risk of injury from bursting hoses during pressure applications due to excessively high temperatures

When operating the pump in pressure applications, hoses that are not designed for the head temperatures of the pump at the respective operating point could become porous and burst.

- Use temperature-resistant pressure hoses at the pneumatic connections.
- Wear protective equipment if necessary (e.g., safety gloves, hearing protection).

**WARNING**

Injury to eyes

Coming too close to the inlet/outlet of the pump may result in injury to the eyes due to the present vacuum/operating pressure.

- Do not look into the pump inlet/outlet during operation.

- Only operate the pump in accordance with the operating parameters and operating conditions described in Chapter 5 *Technical data*.
- Ensure the proper use of the pump (See Chapter 3.1 *Proper use*).
- Eliminate the possibility of improper use of the pump (see Chapter 3.2 *Improper use*).
- Observe the safety notices (see Chapter 4 *Safety*).

Risk of pump head bursting due to excessive pressure increase

- Do not exceed the maximum permissible operating pressure (see 5 *Technical data*).
- Monitor the pressure during operation.
- If the pressure exceeds the maximum permissible operating pressure of the pump: immediately switch off the pump and remedy the fault (see Chapter Troubleshooting).
- Only throttle or regulate the air or gas quantity on the suction line to prevent the maximum permissible operating pressure from being exceeded.
- If the air quantity or gas quantity on the pressure line is throttled or regulated, make sure that the maximum permissible operating pressure at the pump is not exceeded.
- Ensure that the pump outlet is not closed or restricted.

**DANGER**

Risk of dangerous gas mixtures during pump operation

Depending on the medium being transferred, breakage of the media-contacting components can result in a dangerous mixture if the medium mixes with the air in the compressor housing or the surroundings.

- Before using a medium, check the compatibility of the media-contacting components (see 5 *Technical data*) with the medium.

i Excessive pressure, with all of the associated hazards, can be prevented by means of a bypass line with a pressure relief valve between the pressure side and suction side of the pump. Further information is available from KNF Customer Service (contact data: see www.knf.com).

Pump standstill → Establish normal atmospheric pressure in the lines while the pump is at a standstill (relieve pump pneumatically).

9.2 Information on switching the pump on and off

Switching on the pump

- Make certain that the specified pressure or that the specified vacuum (see 5 *Technical data*) is not exceeded in the lines when switching on.
- Switch on the pump by pushing the rotary/push knob (see Fig. 2/3, see Fig. 5).


Fig.5: Switching the pump on/off


Fig.6: Setting the flow rate

Setting the flow rate

The pump speed can be varied with the rotary/push knob. The flow rate can be set in this way (see Fig. 6).

- The speed setting of the rotary/push knob is retained when the pump is switched off.

Switching off/decommissioning pump

- When transferring aggressive media, flush the pump before switching off to extend the service life of the diaphragm (see Chapter 10 *Servicing*).
- Switch off the pump with the rotary/push knob (see Fig. 5).
- If possible, establish normal atmospheric pressure (relieve pump pneumatically to extend the service life of the pump).
- Pull the plug of the power supply of the pump from the grounded socket.

WARNING

Automatic restart after interruption of the power supply

If the power supply is interrupted, the pump automatically restarts.

- Take appropriate measures if necessary.

10 Servicing

ESD-sensitive parts (ESDS)

Failure to observe the ESD protection provision acc. to IEC 61340-5-1 can result in total or partial damage to the pump.

- Maintenance of the pump may only be performed by a qualified person in an ESD-protected area (EPA) acc. to directive IED 61340-5-1.

10.1 Servicing schedule

Risk of explosion from the formation of potentially explosive atmosphere

Leaky connections can result in dangerous potentially explosive atmospheres.

- Ensure that diaphragms and valve plates/seals are installed undamaged, cleanly and correctly.
 - Check the pneumatic connections of the pump for leaks.
 - Work with care during service work.
 - Replace defective parts immediately.
-

Danger of injury through not using genuine spare parts

The functionality of the pump and its safety will be lost, if genuine spare parts are not used.

The validity of the CE conformity is rendered void if genuine parts are not used.

→ Use only genuine spare parts from KNF for servicing work.

Component	Servicing interval
Pump	→ Inspect the pump periodically for external damage or leakage. → Periodically check for noticeable changes to noises and vibrations.
Diaphragm and valve plates/ seals	→ At the latest, replace when the pump flow rate decreases.

Tab.12: Servicing schedule

10.2 Cleaning

During cleaning work, ensure that no fluids enter the interior of the housing.

10.2.1 Flushing the pump

→ Before switching off, flush the pump with air at atmospheric conditions (ambient pressure) for about 5 minutes (if necessary for safety reasons: with an inert gas).

10.2.2 Cleaning the pump

- Only clean the pump with a damp cloth and non-flammable cleaning agents.
- If compressed air is present, blow out the parts.

10.3 Replacing diaphragm and valve plates

- Prerequisites
- Disconnect the pump from mains and ensure that it is voltage-free.
 - Clean the pump and free the pump of hazardous materials.
 - Remove the hoses from the pneumatic pump inlet and pump outlet.
 - If necessary, remove attachment parts such as silencer, suction filter, etc., from the gas connections.

Spare parts/tool

Spare part/tool	Quantity
Spare part set*	1
Torx TX10 screwdriver with torque indicator	1

Tab. 13: Spare parts/tool

*acc. to Chapter 11 Spare parts and accessories

Information on the procedure Diaphragm and valve plates/seals are the only wear parts in the pumps. They are easy to replace.

Valve plates/seals and diaphragm should generally be replaced at the same time. If the diaphragm is not replaced at the same time as the valve plates/seals, the specified output of the pump can no longer be ensured after the maintenance is performed.

**WARNING**

Health hazard due to dangerous substances in the pump

Depending on the medium being transferred, caustic burns or poisoning is possible.

- Wear protective equipment if necessary, e.g., protective gloves, goggles.
- Clean the pump with suitable measures.

**CAUTION**

Risk of burns from hot pump parts

The pump head or motor may still be hot after operation of the pump.

- Allow the pump to cool after operation.

The diaphragm and valve plates/seals are to be replaced in the following order:

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

Fig.8: Loosening housing screws

NOTICE

Make certain that the following steps are performed in an ESD-protected area (EPA).

2. Push the housing cover (2) over the gas connections (see Fig. 9).


Fig.9: Pushing housing cover over gas connections (open pump)

3. Carefully remove the housing cover (2) and place it next to the pump housing (3) (see Fig. 10).
- i** When removing the housing cover (2), make certain that no larger tensile forces are exerted on the cables and that the cables are not damaged


Fig.10: Removing the housing cover

The item numbers within the following work instructions refer to Fig. 11.


Fig. 11: Exploded view

b.) Remove pump head

- Loosen the four head plate screws (4) and remove the head plate (5) together with the intermediate plate (7) from the pump housing.

c.) Change diaphragm

1. Lift the diaphragm (8) on opposing side edges. Then grasp the diaphragm (8) and move the diaphragm (8) to the upper reversal point. Unscrew the diaphragm (8) counterclockwise.
2. Check all parts for soiling and clean the parts if necessary (for further information, see Chapter 10.2 *Cleaning*).

3. Screw the new diaphragm (8) onto the connecting rod (9) (clockwise) and hand tighten the diaphragm (8).

i When screwing in the diaphragm (8), ensure that it is not overtightened.
If the diaphragm is overtightened, there is a risk that it could be damaged.

d.) Change valve plates

1. Separate the head plate (5) from the intermediate plate (7).
2. Remove the valve plates/seals (6) from the head plate (5).
3. Please check the valve seats, intermediate plate (7) and head plate (5) for cleanliness; replace these parts in the event of unevenness, scratches or corrosion (contact your KNF Customer Service in this case).

i Insert the new valve plates/seals (6) in the valve seats of the intermediate plate (7); the valve plates/seals (6) for the pressure and suction side are identical; the same applies for the top and bottom of the valve plates/seals (6).

4. By slightly moving the valve plates/seals (6) horizontally, make certain that they are not under tension.
5. Make certain that the valve plates/seals (6) are centered in the valve seats of the intermediate plate (7).
6. Properly dispose of the replaced diaphragm, valve plates/seals.

e.) Mount pump head

1. Place the intermediate plate (7) with valve plates/seals (6) on the housing.
2. Place head plate (5) on intermediate plate (7) according to the centering.

i Place pump head on the compressor housing according to the alignment of the gas connections.

3. Tighten the screws (4) crosswise (tightening torque: 190-200 Ncm).

f.) Close pump housing

1. Carefully place the housing cover (Fig. 7/2) on the pump housing (3) (see Fig. 12).

i In doing so, make certain that the cables within the pump housing (3) and the circuit board are not damaged.


Fig.12: Mounting the housing cover

2. Push the housing cover (2) over the gas connections until the geometry of the housing cover (2) is aligned with that of the pump housing (3) (see Fig. 13).


Fig.13: Pushing housing cover over gas connections (close pump)

3. Screw in the four housing screws (1) crosswise (tightening torque: 190 – 200 Ncm) (see Fig. 14).


Fig.14: Screwing in housing screws

g.) Final steps

**WARNING**

Risk of injury and poisoning from leaks

- Before recommissioning the pump, check the pump heads and pneumatic connections for leaks. Leaks may lead to poisoning, chemical burns or similar injuries.

-
1. If necessary, screw in attachment parts such as silencer, suction filter, etc., from the gas connections.
 2. Perform a function test.
 - Connect the gas connections to the pump.
 - Connect the pump to the power supply.
 - Test the pump for proper function (among other things, ultimate vacuum).
 - Again disconnect the pump electrically and pneumatically.
 3. Install the pump in the desired installation location according to *8 Set up, installation and connection*.
 - Connect the gas connections to the pump.
 - Connect the pump to the power supply.
 - Test the pump for proper function (among other things, ultimate vacuum).

If you have questions with regard to maintenance, please contact your KNF Customer Service (contact data: see www.knf.com).

11 Spare parts and accessories

i To order spare parts and accessories, please contact your KNF sales partner or KNF Customer Service (contact data: see www.knf.com).

11.1 Spare parts

Spare part set

Spare part set	Order number
Laboport N96	322636

Tab. 14: Spare part set

A spare part set consists of:

Parts	Quantity
Diaphragm	1
Valve plates/ seals	2

Tab. 15: Spare parts

Other spare parts	Order number
Power supply set	323457
Silencer and hose connector set ID6, PP	323456

Tab. 16: Other spare parts

11.2 Accessories

Accessories	Order number
Tripod holder	323484

Tab.17: Accessories


Fig.15: Pump with tripod holder

12 Troubleshooting

Danger to life from electric shock

- All work on the pump may only be performed by an authorized specialist.
- Disconnect the pump power supply before working on the pump.
- Check and ensure that no voltage is present.

→ Check the pump (see following tables).

Pump does not transfer	
Cause	Troubleshooting
No voltage in the electrical mains.	→ Check the room fuse and switch it on if necessary.
Overtemperature protection of the pump has tripped	→ Disconnect pump from electrical mains. → Allow the pump to cool. → Determine the cause of the overheating and rectify.
Connections or lines are blocked.	→ Check the connections and lines. → Remove the blockage.
External valve is closed or filter is plugged.	→ Check external valves and filters.
Condensate has collected in the pump head.	→ Separate the source of the condensate from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas).
Diaphragm or valve plates/seals are worn.	→ Replace the diaphragm and the valve plates/seals (see Chapter 10.3 <i>Replacing diaphragm and valve plates</i>).

Tab. 18: Troubleshooting: Pump does not transfer

Flow rate, pressure or vacuum too low	
The pump is not achieving the flow rate specified in the technical specifications or in the data sheet.	
Cause	Fault remedy
Condensate has collected in the pump head.	<ul style="list-style-type: none"> → Separate the source of the condensate from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas). → If present, open the gas ballast and flush the pump head.
There is overpressure on the pressure side and at the same time vacuum or pressure above atmospheric pressure on the suction side.	<ul style="list-style-type: none"> → Change the pneumatic conditions.
Pneumatic lines or connection parts have insufficient cross section or are throttled.	<ul style="list-style-type: none"> → Disconnect the pump from the system to determine the output values. → Eliminate throttling (e.g., valve) if necessary. → Use lines or connection parts with a larger cross section if necessary.
Leaks occur at connections, lines or pump head.	<ul style="list-style-type: none"> → Ensure the correct seating of the hoses on the hose connectors. → Ensure that the connections are correctly mounted. → Replace the leaky hoses. → Eliminate the leaks.
Connections or lines are completely or partially plugged.	<ul style="list-style-type: none"> → Check the connections and lines. → Remove any parts and particles that are causing plugging.
Head parts are soiled.	<ul style="list-style-type: none"> → Clean the head components.
Diaphragm or valve plates/seals are worn.	<ul style="list-style-type: none"> → Replace the diaphragm and the valve plates/seals (see Chapter 10.3 <i>Replacing diaphragm and valve plates</i>).

Flow rate, pressure or vacuum too low	
The pump is not achieving the flow rate specified in the technical specifications or in the data sheet.	
Cause	Fault remedy
Replaced diaphragm and valve plates/seals.	<ul style="list-style-type: none"> → Check the hoses for leaks. → If necessary, carefully tighten the outer screws of the pressure plate crosswise.
Rotary/push knob is not set to max. speed.	<ul style="list-style-type: none"> → Set the rotary/push knob to max. speed.

Tab. 19: Troubleshooting: Flow rate, pressure or vacuum too low

Fault cannot be rectified

If you are unable to identify any of the specified causes, contact KNF Customer Service or send the pump to KNF Customer Service (contact data: see www.knf.com).

Preparing for return

1. Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with inert gas) to free the pump head of dangerous or aggressive gases (see Chapter 10.2.1 *Flushing the pump*).
2. Remove the pump.
3. Clean the pump (see Chapter 10.2.2 *Cleaning the pump*).
4. Send the pump together with completed Health and Safety Clearance and Decontamination Form to KNF, stating the nature of the transferred medium.

13 Returns

Prerequisite for repairing a pump by KNF is a completed Decontamination Form.

This is made available on the KNF website as a download.

- To find the form, select your country on the overview page (www.knf.com).

You can find the Decontamination Form in the download area.

If you have questions, please contact your sales partner (contact data: see www.knf.com).

14 Index

A	
Accessories	41
C	
Cleaning	29
Commissioning	24
Connecting rod	17
D	
Decontamination Form	44, 45
Diaphragm	17, 35
Dimensions	15
Directives	13
Disposal	12
E	
Eccentric	17
Electrical data	14
Environmental protection	12
ESD	28
F	
Flow rate	14
Setting the flow rate	27
G	
Gas connection	14
H	
Head plate	35
Head plate screw	35
Health and Safety Clearance Form	44
Hose connector	4, 16
Housing cover	33
Housing screw	33
I	
Inlet	16
Inlet valve	17
Installation altitude	15
Installation location	20
Intermediate plate	35
K	
KNF Customer Service	44
L	
Laboport® N96	4
M	
Maintenance	28
Material	14
Medium	9
Mounting	20
O	
Operating pressure	14
Outlet	16
Outlet valve	17
P	
Parameters	
Operating parameters	9, 24
Transport parameters	19
Performance, pneumatic	14
Personnel	11
Pneumatic data	14
Power supply	4

Connecting power supply	21
Plug insert	4, 21
Power, electrical	14
Pressure application	20
Pressure operation	20
Protection class of pump	15
Pump drive	17
Pump housing	34

R

Recycling	12
Repair	45
Replacing wear parts	31

S

Silencer	4, 16, 20
Spare part	40
Spare part set	40
Spare parts	31
Standards	13
Startup parameters	15

T

Temperature	
Ambient temperature	15
Media temperature	15
Operating temperature	15
Storage temperature	19
Transport temperature	19
Tools	31
Transfer chamber	17
Transport	18
Tripod holder	41
Troubleshooting	42

U

Ultimate vacuum	14
-----------------	----

V

Vacuum operation	20
Valve plate	36

W

Warning notice	6
Weight	15

KNF worldwide

You can find our local KNF partners at: www.knf.com

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