



Pump Systems Inc.  
SyringePump.com



Programmable and Customizable Syringe & Peristaltic Pumps

## NE-1000 Single Syringe Pump

**High Pressure Syringe Pump**  
NE-1010

**Continuous Infusion  
Syringe Pump System**  
Dual-NE-1000

**Microfluidics Single  
Syringe Pump**  
NE-1002X



### The NE-1000 Series Syringe Pump Features

Holds 1 syringe from the smallest size available up to 60 mL. A 140 mL syringe can be filled up to 120 mL. NE-1000 & Dual-NE-1000 pumping rates are as low as 0.73  $\mu\text{L/hr}$  with a 1 mL syringe or as high as 35.33 mL/min with a 60 mL syringe. NE-1010 pumping rate are as low as 1.459  $\mu\text{L/hr}$  with a 1 mL syringe or as high as 127.2 mL/min with a 60 mL syringe. NE-1002X pumping rate as low as 0.007 nL/hr with a 0.5  $\mu\text{L}$  syringe or as high as 1448  $\mu\text{L/min}$  with a 60 mL syringe.

### ***Built for Automation***

- Operates stand-alone or from a computer.
- Infuses and withdraws.
- Applications range from simple infusions to complex pumping programs.
- Programmable preset protocols.
- Program up to 41 pumping phases: change pumping rates, set dispensing volumes, insert pauses, control. and respond to external signals, sound the buzzer.
- RS-232 and TTL logic control interfaces.
- Monitor, network, and control up to 100 pumps from one computer.
- Motor stall detection; Non-volatile memory of all parameters and programming.

Connecting two pumps with a dual cable creates a **Dual Pump System** allowing for seamless continuous infusion or emulsification. Upgrade to the X or X2 advanced firmware versions for gradient pumping and increased program memory. Dispensing accuracy of +/-1 %. Unlimited lifetime technical support and Two-year warranty.

***Plus, many, many more features!***

**Not For Clinical Use on Humans**

# NE-1000 Single Syringe Pump Maximum and Minimum Flow Rates

Syringe Manufacturer	Syringe (mL)	Inside Diameter (mm)	Maximum Rate (mL/hr)	Minimum Rate (µL/hr)	Maximum Rate (mL/min)
Air-Tite	0.5	3.8	34.7	0.477	0.578
	50	29.2	2049	28.16	34.15
B-D	1	4.699	53.07	0.73	0.884
	3	8.585	177.1	2.434	2.952
	5	11.99	345.5	4.748	5.758
	10	14.43	500.4	6.876	8.341
	20	19.05	872.2	11.99	14.53
	30	21.59	1120	15.4	18.67
HSW Norm-Ject	1	4.7	53.09	0.73	0.884
	3	9.7	226.1	3.107	3.769
	5	12.5	375.5	5.16	6.259
	10	15.9	607.6	8.349	10.12
	20	20.1	971	13.35	16.18
	30	22.6	1227	16.87	20.45
Monoject	1	5.74	79.18	1.088	1.319
	3	8.941	192.1	2.64	3.202
	6	12.7	387.6	5.326	6.46
	12	15.72	593.9	8.161	9.899
	20	20.12	972.9	13.37	16.21
	35	23.52	1329	18.27	22.15
Terumo	1	4.7	53.09	0.73	0.884
	3	8.95	192.5	2.646	3.208
	5	13	406.1	5.581	6.769
	10	15.8	600	8.244	10
	20	20.15	975.8	13.41	16.26
	30	23.1	1282	17.63	21.37
Poulten & Graf (Glass)	1	6.7	107.8	1.483	1.798
	2	8.91	190.8	2.622	3.18
	3	9.06	197.2	2.711	3.288
	5	11.75	331.8	4.559	5.53
	10	14.67	517.2	7.107	8.62
	20	19.62	925.2	12.72	15.42
Steel Syringes	1	9.538	218.6	3.005	3.644
	3	9.538	218.6	3.005	3.644
	5	12.7	387.6	5.326	6.46
	8	9.538	218.6	3.005	3.644
	20	19.13	879.5	12.09	14.65
	50	28.6	1965	27.01	32.76
	100	34.93	2932	40.29	48.87
	200	44.75	4813	66.13	80.21

Syringe Manufacturer	Syringe (µL)	Inside Diameter (mm)	Maximum Rate (µL/hr)	Minimum Rate (µL/hr)	SGE Syringe (mL)	Inside Diameter (mm)	Maximum Rate (mL/hr)	Minimum Rate (µL/hr)
SGE (Gas Tight)	5	0.343	282.7	0.004	0.25	2.303	12.74	0.176
	10	0.485	565.3	0.008	0.5	3.257	25.49	0.351
	25	0.728	1273	0.018	1	4.606	50.99	0.701
	50	1.03	2549	0.036	2.5	7.284	127.5	1.752
Hamilton Microliter (Glass)	100	1.457	5102	0.071	5	10.3	254.9	3.504
	0.5	0.103	25.49	0.001	10	14.57	510.2	7.01
	1	0.146	51.23	0.001	25	23.03	1274	17.52
	2	0.206	101.9	0.002	50	27.5	1817	24.98
	5	0.326	255.4	0.004	100	34.99	2942	40.43

Not For Clinical Use on Humans

## Specifications

<u>Model</u>	<u>Style</u>	<u>Stall Detection</u>	<u>Number of Syringes</u>	<u>Maximum Syringe Size</u>
NE-1000	Stand-Alone	Yes	1	60 mL; 140 mL partially filled
NE-500	OEM	No	1	60 mL; 140 mL partially filled
NE-501	OEM	No	1	60 mL; 140 mL partially filled

### Mechanical

<b>Motor type:</b>	Step motor
<b>Motor steps per revolution:</b>	400
<b>Motor to drive screw ratio:</b>	15/28
<b>Drive screw pitch:</b>	20 revolutions/”
<b>Micro-stepping:</b>	1/8 to 1/2 depending on motor speed
<b>Advance per step:</b>	0.2126116 $\mu$ m to 0.8504464 $\mu$ m depending on motor speed
<b>Dimensions:</b>	8 3/4” x 5 3/4” x 4 1/2” (LxWxH) (Non-OEM Versions) (22.86 cm x 14.605 cm x 11.43 cm)
<b>Weight:</b>	3.8 lbs. (1.63 kg)
<b>Allen Wrench:</b>	3/32 Hex (Not all models)

### Electrical

<b>Power supply type:</b>	Regulated external wall adapter, power source specific
<b>Power supply output rating:</b>	12V DC @ 1000 mA
<b>Power connector:</b>	2.1 mm, center positive, DC
<b>Voltage at power connector:</b>	12V DC at full load
<b>Amperage:</b>	750 mA at full load

### Operational

<b>Accuracy:</b>	Within 1% error
<b>Reproducibility:</b>	Within 0.1% error
<b>Maximum force:</b>	45 lbs. at minimum speed, 18 lbs. at maximum speed
<b>Syringe inside diameter range:</b>	0.100 to 50.00 mm
<b>Maximum speed:</b>	5.100464828 cm/min
<b>Minimum speed:</b>	0.004204478 cm/hr
<b>Maximum pumping rate:</b>	1699 mL/hr with a B-D 60 mL syringe
<b>Minimum pumping rate:</b>	0.73 $\mu$ L with a B-D 1 mL syringe
<b>Number of Program Phases:</b>	41
<b>RS-232 pump network:</b>	100 pumps maximum
<b>RS-232 selectable baud rates:</b>	300, 1200, 2400, 9600, 19200

### Custom Applications

For specialized and OEM applications contact New Era Pump Systems Inc.  
Custom modifications can be made to the mechanics or the firmware.

**Not For Clinical Use on Humans**