

PH-ABT-NSF-UCBI-0404G-ADA-LH

Product Description

External probe access

These premier built-in undercounter refrigerators are designed in accordance with the NSF/ANSI 456 Standard for Vaccine Storage. With NSF certified, units protect pharmaceuticals at optimal temperatures, preventing waste and allowing for peak delivery.

The ADA compliant, glass door refrigerators utilize microprocessor controllers and feature temperature alarms, remote alarm contacts, LED interior lighting, and probe access ports with included probes. American Biotech Supply Vaccine Storage Refrigerators utilize HFC-free refrigerant for environmental health and energy efficiency.

General Description and Application

Description Single Glass Door Pharmacy/Vaccine Undercounter Refrigerator Built-In ADA Compliant LH

Indoor use only, +18°C to +26°C (+65°F to +78°F), <70% RH Operational environment

Storage capacity 4.6 cu. ft. gross volume

One swing glass door, self-closing, left hinged, non-reversible, magnetic sealed gasket, keyed Door

Shelves Three shelves (two adjustable/one fixed) with guard rail on back

Low profile roller wheels and leveling legs Mounting

Shielded, switched LED lighting, full coverage, balanced spectrum Interior lighting

Rear wall port (1/2") dia.

Forced Air technology, patent pending Airflow management

Insulation Cabinet is foamed-in-place with EPA compliant high density urethane foam

Exterior materials White powder coated steel

Access control Pyxis[®], Omnicell[®] and AcuDose RX[®] compatible

Two (2) years parts and labor warranty, excluding display probe calibration General warranty

Compressor warranty Five (5) years compressor warranty

100 lbs. **Product Weight Shipping Weight** 140 lbs. Rated Amperage 1.74 Amps

NEMA 5-15 plug, 8 to 10 ft typical, conforms to UL471 requirements, Vaccine storage power Power Plug/Power Cord

cord warning label

110-120V AC: 15 A (minimum) Facility Electrical Requirement

Certified in accordance with the NSF/ANSI 456 Standard for Vaccine Storage. UL, C-UL, ETL, C-Agency Listing and Certification ETL listed (either single or dual agency listings) and certified to UL471 standard, hydrocarbon

refrigerant safety.

Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years

certification of calibration, "buffered" probe in the product simulated solution, min/max

memory. F/C switchable, field installable, and visual & audible temp alarm

Pharmacy refrigerator/freezer toolkit and temperature logs

Refrigeration System

Included Accessories

Compressor Hermetic, high performance Refrigerant EPA SNAP compliant, R600a, Isobutane Condenser Hybrid fin and tube with low noise fan

Evaporator Plate wall

Cycle optimized, zero energy Defrost

Performance

Display probe

Simulator ballast

Uniformity¹ (Cabinet air) +/- 0.8°C Stability² (Cabinet air) +/- 1.2°C Maximum temperature variation +/- 1.4°C (Cabinet air)

Temperature rise after an after 8 sec

door openings

Recovery after 3 min door opening Energy consumption

Average heat rejection Noise pressure level (dBA)

Pull down time to nominal operating

Temperature did not exceed 6.4°C at any probe for all required NSF/ANSI 456 testing protocols³

All probes recover to under 8°C within 4.8 min.

1.15 KWh/day⁴

1.57 KWh/day (224 BTU/h)4

43 or less installed

35 min

Controller, Configuration, Alarms and Monitoring

Controller technology Parametric, microprocessor, LED display with 0.1°C resolution

Temperature setpoint range 1°C to 10°C (Setpoint must remain unaltered from the factory setting to remain compliant with

NSF/ANSI 456 performance requirements)

External alarm connection State switching remote alarm contacts

Visual and audible indicators

Calibrated, stainless steel

High / Low temperature, compliant with alarm requirements defined in the NSF/ANSI 456 **Alarms**

Standard for Vaccine Storage

20 ml bottle, glass bead thermal media

Performance data acquired at 22°C ambient, using NSF/ANSI 456 compliant validation ballast probes, empty chamber, during stabilized steady state operation and a DAQ sampling rate of one measurement every 10 seconds

- 1 Uniformity is defined as the maximum variance in temperature across all probes at any point in time over the testing period
- 2 Stability is defined as the maximum variance in temperature experienced by any single probe over the testing period
- 3 Temperature performance for all loaded and unloaded door opening protocols, all alarm, controller and probe requirements as defined in the NSF/ANSI 456 standard for vaccine storage
- 4 Data per Energy Star test results or equivalent testing and calculation. Heat rejection based on daily averages, not continuous operation. Performance exceeds Energy Star requirements.

Product Data Sheet

Undercounter 4.6 cu. ft. Built-in Glass Door Vaccine Refrigerator Left Hinged ADA NSF/ANSI 456 Certified

Certifications










*-one or more of these certifications may apply to this unit.


Temperature Probes							
Probe	Ave	Min	Max				
1	4.6	3.5	5.8				
2	4.9	4.3	5.4				
3	5.0	4.4	5.6				
4	4.6	3.4	5.8				
5	5.0	4.6	5.3				
6	5.3	4.7	5.9				
7	4.8	4.2	5.5				
8	5.1	4.5	5.8				
9	4.8	3.9	5.8				
10	4.8	3.9	5.8				
11	5.5	4.9	6.2				
12	5.1	4.6	5.6				
13	4.9	4.3	5.5				
14	4.9	4.0	5.9				
15	5.5	4.9	6.2				




Temperature Charts









Refrigerator Left Hinged ADA NSF/ANSI 456 Certified



Images





Dimensions								
	Width	Depth	Height	Door Swing	Total open Depth			
Exterior	23 7/8"	26"	31 15/16"	23 1/2"	46"			
Interior	19 1/4"	17 1/2"	22"					

