



designed for scientists



TS 25-3 Cleaning tool

/// Data Sheet

Cleaning tool in combination with the dispersers T 25 digital, T 25 easy clean digital/control

Recommended cleaning process in three steps:

1. Pre-cleaning with suitable solvent
2. Cleaning with detergent
3. Final cleaning and rinsing with demineralized water

Scope of delivery: one tube holder, three tubes with caps.



designed for scientists

Technical Data

Weight [kg]	0.661
Permissible ambient temperature [°C]	5 - 40