



designed for scientists



## RSE XL Stirring bar remover

/// Data Sheet

Universal magnetic stirring bar remover for large magnetic stirring bars.

Includes mechanical coupling for removing the IKAFLON 118 float mounted stirring drive.

Length: 1.04 m

Material: PTFE, VMQ (FDA-compliant)

Note: Before removing stirring bars from large vessels, please remove the I-MAG drive unit.



designed for scientists

#### Technical Data

Material in contact with medium	PTFE+VMQ
Depth [mm]	105
Height [mm]	1040
Weight [kg]	1.6
Permissible ambient temperature [°C]	5 - 40

