



designed for scientists



CA.UTTD Centrifuge adapter

/// Data Sheet

Adapter to provide rigid structural support to the ULTRA-TURRAX® disposable tubes for use in post-grinding centrifuge separation.

Using UTTD tubes directly in the centrifuge rotor is not advised without the centrifuge adapter as this prevents centrifugal forces from destroying the tube membrane.

Quantity: 2 pieces

Compatible with 20 ml and 50 ml tubes



designed for scientists

Materials:

Adapter: PETP

Pin : 304SS

Compatibility:

CA.UTTD with 20 ml UTTD tubes – Fit centrifuge rotors with minimum tube diameter opening \varnothing 29mm

CA.UTTD with 50 ml UTTD tubes – Fit centrifuge rotors with minimum tube diameter opening \varnothing 48mm



designed for scientists

Technical Data

Working temperature max. [°C]	230
Sterilization methods	all methods
Width [mm]	29
Height [mm]	33
Permissible ambient temperature [°C]	5 - 40