



designed for scientists



HA.sp.m.0.5 Micro sparger, 5 µm

/// Data Sheet

The micro sparger (gassing pipe) is used for gassing in the medium. With a pore diameter of 5 µm it is particularly suitable for fine gassing of sensitive cell types (e.g. mammalian cells). This sparger can be used in 0.5 l glass vessels.



designed for scientists

Technical Data

Material in contact with medium	stainless steel 1.4404
Sparger	Micro sparger
Sparger, outer diameter [mm]	8
Sparger, inner diameter [mm]	6
Sparger pore size [μm]	5

