



designed for scientists



TW.IX.28 Vessel holder

/// Data Sheet

Vessel holder for nine sample vessels with diameter of 28 mm.

The holder can be easily placed on the TW.IX multi-position attachment to simplify sample handling and to ensure that the sample vessels do not fall over during mixing.

When mixing at high speeds, it is recommended to use the TW.IX.ANTISLIP mat. This effectively reduces the vibrations that occur when sample vessels come into contact with the vessel holder. Noise generation is thus significantly reduced.



designed for scientists

Technical Data

Hole diameter [mm]	28
Dimensions (W x H x D) [mm]	99 x 47.5 x 101
Weight [kg]	0.2

