



designed for scientists



HA.Is.10 Lid stand

/// Data Sheet

This sturdy lid stand can be used for 5 and 10 l glass vessels. You can conveniently hang the lid on the stand, so that the individual components can be easily mounted or dismantled. In addition, you have the option to put the motor and sensors in designated places, so that the workplace remains tidy. The pressure spring tripod rod also makes it easy to move the tripod up and down.



designed for scientists

Technical Data

Telescope stand stroke [mm]	390
Maximum load [kg]	9.0
Dimensions (W x H x D) [mm]	460 x 1010 x 420
Weight [kg]	5.8
Permissible ambient temperature [°C]	5 - 40

