



designed for scientists



## C 200.2 Conversion kit for C 5012

/// Data Sheet

Customers who purchased the C 200 as a standard package with a C 5010 decomposition vessel require the C 200.2 conversion kit in order to also use the catalytic activated and halogen resistant decomposition vessel C 5012.

The kit contains three feet to raise the shorter decomposition vessel C 5012 up in the calorimeters inner vessel as well as a plate to use with the C 248 oxygen filling station.



designed for scientists

### Technical Data

Weight [kg]	0.069
Permissible ambient temperature [°C]	5 - 40

