

radwag.com

AS 310.5Y Analytical Balance





AS 310.5Y Analytical Balance

The drawings, photos and graphics used are for illustrative purposes only.

Functions

Wi-Fi

Q	Autotest		Dosing	%	Percent Weighing	**	Parts counting
MAX	Peak hold		Formulation	/	Newton unit measurement	<u>.al</u>	Statistics
- <u>0K</u> +	Checkweighing	4	IR sensors	GLP	GLP Procedures	4	Animal weighing
1	Pipettes Calibration	≋	Air density correction	ρ	Density determination		Differential weighing
	Ambient conditions monitoring	SQC	Statistical Quality Control	е	Packaged Goods Control		ALIBI Memory

Datasheet

	AS 310.5Y Analytical Balance			
Metrological parameters				
Maximum capacity [Max]	310 g			
Ainimum load	-			
Readability [d]	0,1 mg			
/erification unit [e]	-			
Fare range	-310 g			
Standard repeatability [5% Max]	0,07 mg			
Standard repeatability [Max]	0,1 mg			
Standard minimum weight (USP)	140 mg			
Standard minimum weight (U=1%, k=2)	14 mg			
Permissible repeatability [5% Max]	0,12 mg			
Permissible repeatability [Max]	0,15 mg			
inearity	±0,3 mg			
Stabilization time	2,5 s			
Adjustment	internal (automatic)			
DIML Class	-			
Physical parameters				
eveling system	semi-automatic - LevelSENSING			
isplay	10" graphic colour touchscreen			
Veighing chamber doors	manual			
Delivery components	Balance, weighing pan, weighing pan shield, bottom cover, power supply.			
Veighing chamber dimensions	190×190×227 mm			
Veighing pan dimensions	ø100 mm			
Packaging dimensions	600×400×550 mm			
Net weight	7,3 kg			
Gross weight	9,3 kg			
Construction				
Protection class	IP 43			
Components and software				
Database capacity	7			
Features of use				
Touch-free operation	2 IR Sensors			
Communication interface				
Communication interface	2×USB-A, USB-C, HDMI, Ethernet, Wi-Fi, Hotspot			
Electrical parameters				
Power supply	Adapter: 100 – 240V AC 50/60Hz 0,6A; 12V DC 1,2A Balance: 12 – 15V DC 0,8A max			
Power consumption max.	4 W			
Environmental conditions				
Operating temperature	+10 ÷ +40 °C			
Ambient conditions monitoring (option)	THBR 2.0 System, THBR BOX, THB P, THB W, THB S			
Relative humidity	40% ÷ 80%			

Repeatability is expressed as a standard deviation from 10 weighing cycles. Stabilization time depends on the ambient conditions and the dynamics of weighing pan loading; specified for FAST profile. 1 Barcode scanners, available as weighing instrument accessory, communicate with the instrument via RS232 interface exclusively.



Accessories

Antivibration Tables
Holders for laboratory flasks
Power Adapters
RS 232, RS 485 cables
Cigarette lighter receptacle power supply cables
Density determination KIT
Additional modules
USB cable (scale - printer)
Professional Weighing Tables
Protective cover for balances
Barcode scanners
Holders for test tubes and filters

Workstation for Pipettes Calibration
THBR 2.0 System - Ambient Conditions Monitoring
Weighing dishes
Antistatic ionizer
Receipt Printer
Fingerprint Reader
RS 232, RS 485 cables
Protective cover for balances
Under-pan weighing
RS 232 cables (scale - printer)
RS 232 – RS 485 Converter

Software

E2R System Label Editor R02 R-LAB RADWAG Development Studio RAD-KEY RADWAG Remote Desktop Scales Editor 2.1

Device dimensions

