

Simport **UriFlex™**

**Free up storage space with the UriFlex™
UriFlex™ 24-Hr Flexible Collection Bag**

The collection bag is flat when empty and easy to store.
When in use, the free-standing design stays upright and is stable, requiring minimal space in a refrigerator.



Simport
Scientific inc.
Since 1975
A family owned & Canadian company
www.simport.com



Simport
Since 1975 *Scientific inc.*

A family owned Canadian company

simport.com

24 HR URINE COLLECTION
 PATIENT INFORMATION
RECOLECCIÓN DE ORINA 24 HORAS
 REGISTRO DE PACIENTE

PATIENT NAME — NOMBRE PACIENTE _____

REFERENCE # — # REFERENCIA _____ DATE OF BIRTH - FECHA DE NACIMIENTO _____

PRESERVATIVE ADDED — CONSERVANTE AÑADIDO _____

TOTAL VOLUME COLLECTED — VOLUMEN FINAL _____

COLLECTION START TIME — FECHA DE INICIO :

HOUR — HORA _____ DATE — FECHA _____

COLLECTION END TIME / FECHA DE FINALIZACIÓN :

HOUR — HORA _____ DATE — FECHA _____

Reduces the amount of waste when compared to rigid containers

Rotating pour spout
 Polypropylene closure and
 Polyethylene spout

Anatomically contoured funnel is
 ideal for people with health
 conditions or impairments

Patient identification area

Large handle

50 mm wide opening

Made of UV resistant materials appropriate for testing light-sensitive analytes

Strong, thick co-extruded laminated plastic film that will not rip, tear or wrinkle even under heavy weight or stress.

Exterior made of PET
 Interior made of PE

Gusset bottom panel allows the bag to stay in self-standing position when full



Dripless pouring with controlled flow

Dimensions: - Empty: 33 cm x 26.5 cm H.
 - Full (approximate) : 29 cm x 16 cm x 25 cm H.

- 50 mm wide
- Tare weight of 64 g (± 1 g)
- No graduations, liquid volume must be determined by weight

Cat. #	Description	Vol. (L)	Qty/Cs
B351-4	UriFlex™ 24-HR flexible collection bag	4	30
B351-4FU	Funnel for UriFlex™ 24-Hr		30



2588 Bernard-Pilon
 Beloeil, Quebec
 J3G 4S5
 Canada

T-450 464-1723
 F-450 464-3394
 info@simport.com
 www.simport.com

Distributed by: