

E

6" CHART RECORDER

JOB NAME: _____

ITEM NO.: _____

NO. REQUIRED: _____

PART NUMBERS:

- 95-6233 (120V)
- 95-6234 (230V)

SIZE: 7-11/32" Wide x 8-43/64" Height x 4-5/16" Deep
(186.5mm x 220mm x 109.5mm)

DESCRIPTION:

This pressure sensitive recorder is specifically made for use with the Market Forge Sterilmatic Sterilizer (*autoclave*):

Domestic Sterilizer Models:

- STM-E / 95-2678
- STM-EL / 95-3441

Export Sterilizer Models:

- STM-EX / 95-2902
- STM-ELX / 95-2903

Its use is for monitoring and recording the sterilizers operating temperature throughout its sterilizing cycles.

One recording chart can be used for a full 24 hour period. Inherent reliability and precision assure accurate results. The recorder can be wall mounted or surface mounted.

ELECTRICAL REQUIREMENTS:



Unit will be rated at 120 or 230VAC, 40mA
120V, 1 Phase, 50/60Hz
230V, 1 Phase, 50/60Hz

STANDARD FEATURES:

- Records in °F (*Fahrenheit*) - 160°F / 280°F
- Records in °C (*Celsius*) - 70°C / 135°C
- Box of 50 Recording Charts.
- Pressure Sensitive Stylus (*inkless*).
- 6' Cord and Plug
- Temperature Sensor – RTD 1/4"OD x 8" 300 Series Stainless Steel, 1/4" x 1/2" Stainless Steel compression fitting, and 7' cable.
- UL Approved.

EXTRA CHARTS:*(OPTIONAL AT EXTRA COST)*

- Box of 50, 6" pressure sensitive charts, range from 160°F to 280°F.
- Box of 50, 6" pressure sensitive charts, range from 70°C to 135°C.



The manufacturer reserves the right to modify materials and specifications without notice. It is our policy to build equipment which is design certified by U.L. However, a continuing program of product improvement makes it necessary to submit new models to the agencies as they are developed and consequently not all models bear the appropriate agency labels at all times.

SPEC SHEET: S-6048A
08/12LABORATORY EQUIPMENT
6" Chart Recorder

35 Garvey Street • Everett • MA • 02149-4403
Tel: (617) 387-4100 • Toll Free: (866) 698-3188
Fax: (617) 387-4456 (MA and Overseas) • Fax: (800) 227-2659 (Ex. MA)
E-mail: custserv@mfi.com • www.mfi.com