


# Nanodyne Replacement Illuminator for Leica M655 Surgical Microscope - Included Items

## Additional Items Included But Not Shown:

PN 10337 Hex Key 2mm short arm (to attach Nanodyne illuminator)

PN 11709 Hex Key 1.5mm long arm (for M2 latch screws)

PN 11115 SHCS M2 x 0.4 x 4mm (2 pcs) (for attaching spring latch to new illuminator base)



PN 11661 Leica M655 illuminator Assy

PN 10362 M3 x 0.5 x 6mm  
Button Head Cap Screw  
(4 pcs)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11663 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.



PN 11000 Power Cord Extension 1.35mm ID x 3.5mm OD RA plug to jack 6 foot. (use if needed for more length)

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	3/18/2020
PN 11708 Leica M655 illuminator System	REV 1
SHEET 1	OF 1