



Snap-Cap Microcentrifuge Tubes



Various sizes of the Snap-Cap Microcentrifuge Tubes

The **Snap-Cap Microcentrifuge Tube** gives a firm, dependable seal, and its tough virgin-polypropylene construction withstands the rugged G-forces of high-speed centrifugation. Use it with confidence for specimen collection, mixing, electrophoresis sampling/handling, enzyme assays, radioassays, immunoassays and many other applications.

The **Snap-Cap Microcentrifuge Tube** is available in five sizes, 250µl, 400µl, 500µl, 1.5ml and 1.8ml. All sizes have attached, tight-fitting Snap-Caps for optimal sealing during sample processing and centrifugation. Colors: Natural (white) for 250µl, 400µl, and 500µl. The 1.5 ml and 1.8ml sizes are available in natural and four additional colors for color-coding.

The **Snap-Cap Microcentrifuge Tube's** smooth, conical interior lets you remove your entire sample with ease, and without material hangup. You can even cut it to remove a spun pellet. The Snap-Cap is attached to each tube, and the entire unit is autoclavable at 250°F at 15 psi for 15 minutes. Fits Beckman™, Eppendorf™ and all other microcentrifuges.

Ordering Information

Catalog #	Description	Pkg. Qty	Pkgs./ Case
4003	Micro-Centrifuge Tube 400µl, Polypropylene , White	1000	10
4004	Micro-Centrifuge Tube 400µl, Polyethylene, White	1000	10
4005	Micro-Centrifuge Tube 500µl, Polypropylene, White	1000	10
4025	Micro-Centrifuge Tube 250µl, Polyethylene, White	1000	10
4150	Micro-Centrifuge Tube 1.5ml, Polypropylene , White	1000	10
4151	Micro-Centrifuge Tube 1.5ml, Polypropylene , Red	1000	10
4152	Micro-Centrifuge Tube 1.5ml, Polypropylene , Yellow	1000	10
4153	Micro-Centrifuge Tube 1.5ml, Polypropylene , Blue	1000	10
4154	Micro-Centrifuge Tube 1.5ml, Polypropylene , Green	1000	10
4155	Micro-Centrifuge Tube 1.8ml, Polypropylene , White	1000	10
4156	Micro-Centrifuge Tube 1.8ml, Polypropylene , Red	1000	10
4157	Micro-Centrifuge Tube 1.8ml, Polypropylene , Yellow	1000	10
4158	Micro-Centrifuge Tube 1.8ml, Polypropylene , Blue	1000	10
4159	Micro-Centrifuge Tube 1.8ml, Polypropylene , Green	1000	10
4000	Micro Tube Tool	4	25

Autoclavable
 at
 250°F
 at
 15 psi
 for
 15 mins

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4340 Redwood Hwy, Suite A1, San Rafael, CA 94903 USA Tel: (415) 472-3777 Fax: (415) 472-3758 info@bioplas.com
 www.bioplas.com