




- **Autofil[®] 2** is available as a complete system or funnel only
- Funnel is compatible with GL45 thread bottles
- **Autofil[®] 2** units include adapter for connecting to vacuum hose
- **Autofil[®] 2** design includes the **NEW EZlabpure[™]** ergonomic bottle shape with molded finger grip
- **Autofil[®] 2** patented GL45 caps provides a convenient gripping surface with a frosted area for labeling
- Made with USP Class VI materials
- Manufactured in a Class 100,000 clean-room
- Manufactured in an ISO 13485 Facility in North America
- Optional, patented docking hardware sold separately




The patented, next generation **Autofil[®] 2** features the **NEW High Flow Foxx Velocity 0.1µm, 0.2µm, and 0.45µm PES membranes**. Perfect for filtering FBS and other viscous solutions. Our membrane is low protein binding and surfactant-free. **Autofil[®] 2** is available in 250mL, 500mL, and 1000mL complete units with the **NEW EZlabpure[™]** media bottle or funnel only bottle-top filters. All units are sterile and manufactured from USP Class VI polystyrene in an ISO 13485 clean room, FDA regulated facility.

Autofil[®] 2 is compatible with the patented **Autofil[®] Pedestal** and **NEW Cradle Ring 2**. These optional docking stations allow for hands free filtering and provide a stationary vacuum hose connection, allowing for efficient set-up and quick use while preventing spills and tissue culture hood contamination.

- 1 Patented single-piece design made entirely from polystyrene material
- 2 Features **NEW High Flow Foxx Velocity 0.2µm PES membrane**
- 3 Forward-facing, large, raised graduation marks for easy reading
- 4 Easy to read label makes identifying membrane type and size easy when removing from packaging
- 5 **EZlabpure[™]** PS sterile media bottle with GL45 thread type



ISO 13485



	Item Number	Volume	Membrane Pore Size	Filter Diameter	Description
Full Assembly	116-71132-FLS	250mL	0.1µm PES	66mm	250mL Full Assembly 0.1µm PES, 12/pk
	116-81132-FLS	500mL	0.1µm PES	90mm	500mL Full Assembly 0.1µm PES, 12/pk
	116-91132-FLS	1000mL	0.1µm PES	90mm	1000mL Full Assembly 0.1µm PES, 12/pk
	11012-FLS	250mL	0.2µm PES	66mm	250mL Full Assembly 0.2µm PES, 12/pk
	11022-FLS	500mL	0.2µm PES	90mm	500mL Full Assembly 0.2µm PES, 12/pk
	11032-FLS	1000mL	0.2µm PES	90mm	1000mL Full Assembly 0.2µm PES, 12/pk
	11412-FLS	250mL	0.45µm PES	66mm	250mL Full Assembly 0.45µm PES, 12/pk
	11422-FLS	500mL	0.45µm PES	90mm	500mL Full Assembly 0.45µm PES, 12/pk
	116-3233-RLS	1000mL	0.45µm PES	90mm	1000mL Full Assembly 0.45µm PES, 12/pk
Funnel Only	116-11132-FLS	250mL	0.1µm PES	66mm	250mL Funnel Only 0.1µm PES, 12/pk
	116-21132-FLS	500mL	0.1µm PES	90mm	500mL Funnel Only 0.1µm PES, 12/pk
	116-31132-FLS	1000mL	0.1µm PES	90mm	1000mL Funnel Only 0.1µm PES, 12/pk
	11512-FLS	250mL	0.2µm PES	66mm	250mL Funnel Only 0.2µm PES, 12/pk
	11522-FLS	500mL	0.2µm PES	90mm	500mL Funnel Only 0.2µm PES, 12/pk
	11532-FLS	1000mL	0.2µm PES	90mm	1000mL Funnel Only 0.2µm PES, 12/pk
	11612-FLS	250mL	0.45µm PES	66mm	250mL Funnel Only 0.45µm PES, 12/pk
	11622-FLS	500mL	0.45µm PES	90mm	500mL Funnel Only 0.45µm PES, 12/pk
	11632-FLS	1000mL	0.45µm PES	90mm	1000mL Funnel Only 0.45µm PES, 12/pk
Optional Accessories	11712-FLS	250mL	-	-	250mL EZlabpure™ Media Bottle, PS, 12/pk
	11722-FLS	500mL	-	-	500mL EZlabpure™ Media Bottle, PS, 12/pk
	11732-FLS	1000mL	-	-	1000mL EZlabpure™ Media Bottle, PS, 12/pk
	1124-FLS	-	-	-	Autofil® Quick-Connect Cradle Ring 2, 1/pk
	1116-RLS	-	-	-	Autofil® Quick-Connect Pedestal, 1/pk



Autofil® is a registered trademark of Foxx Life Sciences.

