


ENDURO GDS II GEL DOCUMENTATION SYSTEM GDS II & GDST II

- **Designed to be User Friendly**
- **No Focusing**
- **Auto-exposure works even with faint bands**

Enduro GDS II is the next generation Gel Documentation System from Labnet. The GDS II comes in either a 302nm or a 365nm version and runs on a Windows based PC (not included). The GDST II comes in either a 302nm or 365nm wavelength version with 470nm epi-blue lights for use with “Safe” Dyes that eliminate need for ethidium bromide. Both versions now come with a storage location on its top to help maximize bench space and provide a place for accessories so you know where to look. Ask your local Labnet distributor for a demo today.

	GDS II	GDST II
Image Resolution	3.2 MP	3.2 MP
Field of View	15 x 20 cm	15 x 20 cm
Standard Illumination	UV (302 or 365nm), epi white	UV (302 or 365nm), epi-white, epi-blue
System Controls	Windows 10 PC or laptop (which is not included)	Windows 10 Tablet (built-on)
USB 2.0 ports	2	2
Dimensions	12 x 15 inches (31 x 38 cm)	12 x 15 inches (31 x 38 cm)
Certifications	cTUVus; CE	cTUVus; CE
Computer Specifications	PC or Laptop not included but must meet following minimum specifications: <ul style="list-style-type: none"> ● Windows 10 ● Intel® processor ● 1.4 GHz processor speed ● 2 GB RAM ● 16 GB free hard disk space ● 2 USB 2.0 or higher ports 	Tablet Included

GDS II & GDST II



470nm Epi Blue Lights
Ideal for DNA "Safe Dyes" (GDST II only)



Maximize Valuable Lab Space
with Small Footprint and Storage on Top

Part Numbers

Cat. no.	Description
GDS2-1302	Enduro GDS II with 302nm UV
GDS2-1365	Enduro GDS II with 365nm UV
GDST2-1302	Enduro GDS Touch II with 302nm UV
GDS-TL1D-1	Total Lab 1D single user licence for quantitative analysis
GDS2-GCS	Gel Cutting Surface
GDS2-GCT	Gel Carrying Tray
GDS2-OG	Orange Goggles
GDS2-WLCS	White light Conversion Screen
GDS2-USBS	USB Splitter
GDS2-USBMK	USB Mouse/Keyboard
GDS-12021	Mitsubishi P-95DW Thermal Printer
GDS-12022	High Gloss Thermal Paper
GDS2-BLCS	Blue Light Conversion Screen

APPLICATIONS:

- Visible protein stains like Coomassie blue or silver stain
- Designed for DNA
 - Ethidium bromide
 - Safe Dyes such as GelGreen® and SYBR™ safe

FEATURES

- region of interest function allows optimal exposure in less than 30 seconds. No more multiple exposures to get good images.
- Transilluminator good for 30,000 hours so that no longer need to buy replacement UV bulbs
- 16bit Tiff files ideal for publication
- Included fully featured qualitative software ideal for creating publication and notebook quality images